

# Pen comparison

	Prefilled pens		Durable pens		
Pens	FlexTouch®	FlexPen®	NovoPen® 5	NovoPen® 4	NovoPen Echo®
Minimum dose (units)	1 (100 U/ml pen) 2 (200 U/ml pen)	1	1	1	½
Maximum dose (units)	80 (100 U/ml pen) 160 (200 U/ml pen)	60	60	60	30
Colour differentiation	✓	✓	Comes in two colours and a variety of skins that can be applied to the pen	Comes in two colours and a variety of skins	Comes in two colours and a variety of skins that can be applied to the pen
Memory function			✓		✓
No dose button extension	✓				
Easy-to-read dose display	✓	✓	✓	✓	✓
Dose dial clicks forward and backward	✓	✓	✓	✓	✓
End-of-dose click	✓		✓	✓	✓
Half unit dosing					✓
Compatible with most needles	✓	✓	✓	✓	✓

# Needle comparison

Needles	NovoFine®	NovoFine® Plus	NovoTwist®	NovoFine® Autocover®
<b>Needle size</b> (G* and mm)	32G 4 mm 32G 6 mm 31G 6 mm 30G 8 mm	32G 4 mm	32G 5 mm 30G 8 mm	30G 8 mm
<b>Bayonet fitting</b> (just-twist attachment)			✓	
<b>Screw-on attachment</b>	✓	✓		✓
<b>Silicone coated</b> for injection comfort	✓ <small>(32G electro-polished)</small>	✓ <small>(32G electro-polished)</small>	✓ <small>(32G electro-polished)</small>	✓
<b>ETW or Superflow™ technology</b> makes your injection quicker and easier†	✓	✓	✓	
<b>Integrated glue tower</b> reduces the risk of needle breakage	✓	✓	✓	
<b>Compatible with all Novo Nordisk diabetes pens</b>	✓	✓	✓	✓
<b>Compatible with most major diabetes pens</b>	✓	✓		✓
<b>Concealed needle</b>				✓
<b>Safety shield</b> locks in place after use				✓

\*The higher the number, the thinner the needle.

†The thinner needle wall allows for a larger inner diameter, reducing the effort needed to inject your insulin.